

IMPORTANT: The Humphrey ADE-circle System - Fresh Gas Flows (FGF) Mar 2018

Anaequip-Vet UK, 2 Millstream Bank, Worthen, Shrewsbury, SY5 9EY Tel 01743-891140 email: sales@anaequip.com.

ONLINE SHOP www.anaequip.com for all ADE spares, tubes etc. **Also for education material.**

General comments: (for more copies, in-depth detail, set up, education & future webinars → www.anaequip.com)

- **The ADE is significantly more efficient than all other systems with its new 4-phase exhaust valve. Design changes reduce fresh gas usage by up to 80%; but dialled vapour % remains unaltered.**
- For most healthy animals 50kg or less, a maintenance flow of just 500 ml/min will be enough!
- **The minimum oxygen flow recommended is 300 ml/min.** This allows for small leaks & vaporizer inaccuracy.
- Read "Routine Maintenance Instructions" supplied for **leak tests**, soda lime refills & general info.
- If using **nitrous oxide**, **special conditions apply especially when recycling through soda lime.**
- **Lever position: UP** for **spontaneous & manual IPPV ventilation.** Place tip of finger in exhaust valve screw cap to hold it closed and squeeze bag. **Ventilator mode** – lift/shift lever **DOWN** (details on website).
- **A premedication with analgesic drug** (eg buprenorphine) greatly assists a smooth anaesthetic.
- If a side stream **capnograph or vapour monitor** that samples gas is being used, the **sample flow** (often 100-200 ml/min) **must be added** to the fresh gas flow recommendations below.

Basic Recommendations: All cats, dogs & other small animals up to 7 kg

Circuit factor = 0.5 – 0.7: FGF 70-100 ml/kg/min. The ADE system really is highly efficient!

- **If using nitrous oxide, use as normal: the O₂ flow must be at least 30% of the mixture ie ratio of 2N₂O:1O₂**
- Remove canister - No recycling through soda lime
 - use ½ litre reservoir bag.
- **Leak test** is essential before **every anaesthetic!**
 - Exhaust valve: **fully open** after leak test.
- Connect to patient – check ET tube doesn't leak: **Fill system for 1 minute at 2 l/min with oxygen and vapour** to wash out non-anaesthetic gases in bag and tubing. This considerably speeds up induction.
- Continue with induction flows and reduce to maintenance levels on approaching level of anaesthesia required.
- Adjust vaporizer % on **clinical signs (all animals vary)**, but **keep flows as recommended below.**
- For **quick** adjustments in vapour %, empty bag, refill with new % at high flows & **reduce** when bag full.

| Weight Kg | Induction fresh gas flow pre-fill system with O ₂ + vapour (set at twice % for maintenance) | Maintenance flow premed with analgesic healthy animal | Maintenance flow No analgesia in premed or compromised animal |
|-----------|--|---|---|
| 0-3 kg * | 500 ml/min | 300 ml/min (minimum) | 500 ml/min |
| 3-4 | 500 | 300 | 500 |
| 4-5 | 500 | 350 | 500 |
| 5-6 | 600 | 400 | 600 |
| 6-7 | 700 | 500 | 700 |

* For birds & very small animals below 1kg, connect ADE system to a conical mask or induction box.

Basic Recommendations: Animals over 7 kg up to 100kg

Circuit factor = 0.1 – 0.3: FGF 10 - 30 ml/kg/min. Again the ADE-circle really is highly efficient!

- **If using nitrous oxide, use a 50/50 O₂/N₂O mixture. Set minimum flow at 300ml each = 600 ml/min**
- Attach canister – change soda lime if exhausted
 - use 1 litre reservoir bag (2 litre over 50kg)
- **Leak test** is essential before **every anaesthetic!**
 - Exhaust valve: **fully open** after leak test.
- Connect to patient – check ET tube doesn't leak: **Fill the system for 2 minutes at 3 l/min with oxygen and vapour** to wash out the 3 litres of non-anaesthetic gases in the bag and canister. **This considerably speeds up induction.**
- Continue with induction flows and reduce to maintenance levels on approaching level of anaesthesia required.
- Adjust vaporizer % on **clinical signs (heavy & fat animals need more)**, but **keep flows as recommended.**
- For **quick** adjustments in vapour %, empty bag, refill with new % at high flows & **reduce** when bag full.

| Weight Kg | Induction flow: 30 ml/kg/min. Set vapour % at ± 1.5 times maintenance % to ensure the system is filled with O ₂ + vapour | Maintenance flow healthy animal Reduce flow to 10 ml/kg/min Adjust vapour % as required | Maintenance flow compromised animal use higher flows at 30 ml/kg/min |
|--------------|---|---|--|
| 7-15 kg | 500 ml/min | 300 ml/min (minimum) | 500 ml/min |
| 15-20 | 600 | 300 | 600 |
| 20-30 | 900 | 300 | 900 |
| 30-40 | 1200 | 400 | 1200 |
| 40-50 | 1500 | 500 | 1500 |
| 50-60 | 1800 | 600 | 1800 |
| 60-70 | 2100 | 700 | 2100 |
| 70-80 | 2400 | 800 | 2400 |
| 80-90 | 2700 | 900 | 2700 |
| 90-100 | 3000 | 1000 | 3000 |
| Up to 250 kg | 30 ml/kg/min – use 22mm tubing | Reduce to 10 ml/kg/min | use 3 litre reservoir bag |